

## LCS 708: SLIMMER-8 CHANNEL

LCS 708 Slimmer is the body treatment appliance available in lightweight with all standard accessories for slimming, refirming and anti-cellulite treatments.



### Technical Specifications: -

Operating Voltage	200V ~ 230V, 50Hz
Weight	1.8Kg (Approx)
Room Temperature	10°C ~ 40°C
Moisture	10% ~ 80%
Therapy Mode	Continuous, Intermittent & Pulse
Output Eight Channel	(16 Pads)
Timer	0 to 99 Minutes Digital
Contraction	Adjustable
Relaxation	Adjustable

### Standard Accessories: -

Electrodes (Size: 3" Dia)	16No.
Output Cable	8No.
Velcro Straps	8No.
Mains Cable	1No.
Jelly Bottle	1No.

## **LCS 712: SLIMMER-12 CHANNEL**

**LCS 712 Slimmer** is the body treatment appliance available in lightweight with all standard accessories for slimming, refirming and anti-cellulite treatments.



### **Technical Specification: -**

Operating Voltage	200V~230V, 50Hz
Weight	1 Kg (Approx)
Room Temperature	10°C ~ 40°C
Moisture	10% ~ 80%
Therapy Mode	Continuous, Intermittent & Pulse
Output	Twelve Channel (24 Pads)
Timer	0 To 99 Minutes Digital
Contraction	Adjustable
Relaxation	Adjustable

### **Standard Accessories: -**

Electrodes (Size - 3" Dia)	24No.
Output Cable	12No.
Velcro Straps	12No.
Mains Cable	1No.
Jelly Bottle	1No.

## LCS 716: SLIMMER-16 CHANNEL

LCS 716 Slimmer is the body treatment appliance available in lightweight with all standard accessories for slimming, refirming and anti-cellulite treatments.



### Technical Specification: -

Operating Voltage	200V~230V, 50Hz
Weight	1 Kg (Approx)
Room Temperature	10°C ~ 40°C
Moisture	10% ~ 80%
Therapy Mode	Continuous, Intermittent & Pulse
Output	Sixteen Channel (32 Pads)
Timer	0 To 99 Minutes Digital
Contraction	Adjustable
Relaxation	Adjustable

### Standard Accessories: -

Electrodes (Size - 3" Dia)	32No.
Output Cable	16No.
Velcro Straps	16No.
Mains Cable	1No.
Jelly Bottle	1No.

## LCS 775: VACUUM THERAPY UNIT

**Vacuum Therapy Unit** is very compact & fitted with durable pump; this vacuum therapy is used in slimming centres on abdomen. Hips and thigh areas before starting electro-therapy stimulation.



### Technical Specification:

Operating Voltage	220V, 50Hz
Room Temperature	10°C ~ 40°C
Moisture	10% ~ 80%
Timer	Digital
Display	Analog Vacuum Display

### Standard Accessories: -

Vacuum Cup (Borosil Glass)	2No.
Flexible PVC Tubes	1No.
Mains Cable	1No.

## **LCS 777: MASSAGER (G-5):**



### **OUTSTANDING CHARACTERISTICS: -**

- ☞ This massager combines a vertical motion with an irregular rotary motion, which raises massaging effects.
- ☞ Powerful, irregular vibration gives the same effect as your skillful hands do, enabling the unit to do the job with an almost explosive speed.
- ☞ A massage with heat on is more effective for treatment of stiffness and fatigue. The temperature of heating part raises up to 55°C (131°F). Hot and cool of heating part is changeable with a switch of the unit.
- ☞ Designed for a long use, it can be used with absolute safety, anytime and anywhere.

### **TECHNICAL SPECIFICATION:**

<b>Operating Voltage</b>	<b>AC 220V, 50/60Hz 20W (Hot), 18W (Cold)</b>
<b>Maximum Temperature</b>	<b>55° C (131° F)</b>
<b>Net Weight</b>	<b>1.1 Kg.</b>
<b>Size</b>	<b>300 x 110 x 120 mm</b>
<b>Speed</b>	<b>4500 rpm</b>
<b>Accessories</b>	<b>4 Kinds</b>