

ADVANCE

# AYU-42AB ELECTRO SURGICAL UNIT (400W)

- Light Weight
- Elegant Design
- Simple Operation
- The Reliable Economical

[www.lcsindia.com](http://www.lcsindia.com)



## High Frequency Electro Surgical Unit With Separate Boiler Footswitch

microprocessor control, innovative design, together with high reliability makes 400W High frequency Surgical Units the most favorable choice for doctors.

In the used surgical diathermy, the technical ability of the unit combined with plain safety are of prime importance-both for the patient and the user; There must be no possibility of accidental, undesirable results in the use of Ayunix Electro Surgical Units been developed.

The Ayu 400W-ANL are newly developed Electro surgical Units of particular safety and reliability. The Ayu 400W as the most refined unit of its kind in the world.

### SAFETY IN USE


The Technical conceptions of the Ayu 400W-ANL meet the requirements of highly refined, Electro Surgical techniques all national and international safety standards are conformed to or exceeded. It is not possible for the unit to deliver a higher output than has been set on the control. The micro function allows for an instant reduction in the output for more delicate application such as neurological or micro surgery.

### OUTPUT CHARACTERISTICS

|                    | MODE       | RETED LOAD<br>(OHMS) | MAXIMUM POWER<br>(WATTS) |
|--------------------|------------|----------------------|--------------------------|
| Monopolar<br>Cut   | PURE       | 500                  | 400                      |
|                    | BLEND-I    | 300                  | 275                      |
| Monopolar<br>Coag. | DESSICCATE | 300                  | 250                      |
|                    | FULGURATE  | 300                  | 150                      |
| Bipolar            | CUT        | 100                  | 100                      |
|                    | COAG       | 100                  | 70                       |

Physical Characteristics - L x W x H - 430x405 x 280 mm  
Net Weight - 5Kg 500Grm. Gross Weight - 10Kg. 900Grm.

### Standard Accessory



Patient Plate

Hand Switch

Unipolar  
Handle

Knife Electrode  
(Angled)

Knife Electrode  
(Straight)

Needle Electrode

Ball Electrode  
6mm

Loop Electrode  
5mm

Bipolar  
Forceps

Unipolar & Bipolar  
foot switch



### LIFE CARE SYSTEMS

K-708, Sector-23, Sanjay Nagar,  
Ghaziabad-201002 Delhi NCR  
Cell : +91 9312247417  
Cell : +91 9654875195  
E-mail : sales@lcsindia.com  
Website : www.lcsindia.com